

Scientific Background Information on BrightMEM Allograft

Significance of Research Question/Purpose:

Limbal stem cell deficiency (LSCD) is a devastating disease that accounts for an estimated 15-20% of corneal blindness worldwide.¹ This disease results from loss or dysfunction of LSCs, a population of pluripotent cells that continuously regenerates the transparent epithelium of the cornea throughout life.²⁻⁴ LSCs are found on the ocular surface in the limbus, where they divide and differentiate into transient amplifying cells (TAC), which then migrate centripetally across the cornea and further differentiate into mature corneal epithelium.⁵ Loss of these cells from chemical burns, autoimmune diseases, various congenital disorders, ocular surgeries, and other surface insults can result in partial or total LSCD and associated vision loss.²

Without a healthy population of LSCs to regenerate the corneal epithelium, LSCD patients are at risk for recurrent erosions, persistent epithelial defects (PED), corneal conjunctivalization, corneal scarring, and corneal melting. Standard corneal transplants are ineffective. Instead, transplantation of donor limbal tissue containing donor LSCs (keratolimbal allografts or KLAL) is often necessary. However, because systemic immunosuppression with prednisone, mycophenolate, and tacrolimus is required to prevent rejection of the highly antigenic limbal grafts, KLAL is almost always reserved for only the most severe cases.⁶⁻⁸ For many patients, the risks associated with immunosuppression outweigh the benefits of treatment. Artificial corneas made from plastic polymers such as the Boston keratoprosthesis (KPro) have also been used to treat LSCD because their ability to remain clear in the absence of LSCs and healthy corneal epithelium;⁹ however, the significant rates of glaucoma, corneal melt, implant extrusion, infectious keratitis, and endophthalmitis associated with KPros have also limited their use in mild to moderate disease.¹⁰ As such, current treatment options for LSCD remain limited.

Persistent epithelial defects (PED) are a well-known complication of LSCD that results from the failure of the corneal epithelium to regenerate an injury within 10-14 days.¹¹ PEDs can also occur due to neurotrophic keratopathy, infectious keratitis, chronic exposure keratopathy, or other etiologies where proliferation, migration, and adhesion of the corneal epithelium are impaired. Both LSCD and PEDs are challenging disorders of the corneal epithelial regeneration where therapies to promote durable corneal re-epithelialization remain lacking.

The limbal niche is the microenvironment surrounding the LSCs and limbal TACs that is critical for maintaining their survival and proliferative potential under physiologic conditions.^{12,13} The limbal niche is extremely complex and includes the limbal basement membrane, the limbal stroma, the microstructures of the limbus (such as the limbal crypts), as well as the limbal nerve plexus, microvascular, melanocytes, dendritic cells, and mesenchymal stem cells.¹⁴⁻¹⁹ Extracellular signals from this microenvironment are critical to the normal function and maintenance of pluripotent stem cells.^{20,21}

Without the support of a niche-like environment, LSCs are often dysfunctional or have limited survival over time. This is evidenced by the fact that, in contrast to KLAL, where limbal tissue is transplanted with LSCs, transplantation of cultured LSCs alone often results in poor long-term survival of the transplanted LSCs. Even in the absence of rejection, transplanted cultured LSCs have been reported to have limited survival beyond 9 months.^{22,23}

Similar evidence is also seen in patients with congenital aniridia. In congenital aniridia, limbal niche dysfunction is believed to contribute in part to the later development of LSCD. Due to a PAX6 mutation, patients with congenital aniridia have malformation of their limbus.²⁴ Consequently, although aniridics do not have LSCD early in life, without a normal limbal niche, these patients gradually develop attrition of their LSCs over time and typically present with LSCD in the second or third decades of life.²⁴ As such, the limbal niche is critical for long-term LSC and limbal TAC survival and function.

At the same time, extracellular signals from the limbal niche can be powerful drivers of LSC, TAC, and even corneal epithelial proliferation. In fact, treatments that leverage components of the limbal niche to stimulate epithelial cell proliferation on the corneal surface have been explored extensively. Mesenchymal stromal cells are an important component of the limbal niche. Topical application of mesenchymal stromal cells and mesenchymal stromal cell-derived exosomes have been shown to stimulate corneal epithelial regeneration in both LSCD and PEDs.^{25,26} By activating cell signaling pathways normally stimulated by the extracellular matrix and cells of the limbal niche, limbal and epithelial cell proliferation can be stimulated even outside of the native limbus.

Isolated components of the limbal niche thus holds significant promise as potential therapies for ocular surface diseases. Descemet's membrane (DM), an acellular, naturally occurring basement membrane found on the posterior surface of the cornea, is a promising substitute for limbal basement membrane. DM is routinely isolated and transplanted intraocularly for treatment of diseases such as Fuchs' dystrophy and corneal bullous keratopathy. However, its application on the ocular surface has not been explored.²⁷ DM is optically clear and highly resistant to collagenase digestion.^{28,29} This is evidenced by the fact that Descemetocelles can remain intact for months in the setting of a sterile corneal melt. Taken together, this makes DM attractive as a potential long-term corneal on-lay and substrate for supporting corneal epithelial regeneration.

The anterior fetal banded layer of DM shares key compositional similarities with limbal basement membrane, which is a major component of the limbal niche. These similarities include limbus-specific extracellular matrix proteins such as collagen IV that is restricted to the $\alpha 1$ and $\alpha 2$ subtypes, vitronectin, and BM40/SPARC.³⁰⁻³² Of these, vitronectin and BM40/SPARC are known to promote proliferation of LSCs and induced pluripotent stem cells (iPSC) in culture.^{33,34}

Because of this, DM is a promising biological membrane for establishing a pro-proliferative substrate on the corneal surface in patients with LSCD and PEDs.

	Epithelial BM		DM
	Central	Limbus	Stromal face
Type IV collagen chains	$\alpha 3-\alpha 6$	$\alpha 1, \alpha 2, \alpha 5, \alpha 6$	$\alpha 1, \alpha 2$
Laminin chains	$\alpha 1, \alpha 3-$ $\alpha 5, \beta 1, \beta 3, \gamma 1-3$	$\alpha 2-\alpha 5, \beta 1-$ $3, \gamma 1-3$	none
Nidogen-1/entactin-1	+	+	-
Nidogen-2/entactin-2	+	+	-
Perlecan	+	+	-
Netrin-4/ β -netrin	-	-	-
Matrilin-4	+	+	-
Tenascin-C	-	+	-
Fibrillin-1	-	+	-
Type VIII collagen			+
Type XII collagen	+ (both forms)	+ (short form)	-
Type XVIII collagen	+	+	-
Thrombospondin-1	+	-	-
SPARC/BM-40/osteonectin	-*	+*	\pm
Fibronectin	+	+	+
Vitronectin	-	+	+
Endostatin	+	+	
Versican	-	+	

\pm weak staining in some cases.

**Slotzer-Shrehardt, et al. (Exp Eye Res 2007) reported (-) in the central epithelial BM, and (++) in the limbal BM. However, Kabasova, et al. (Invest Ophthalmol Vis Sci 2008) reported weak staining in some cases for epithelial BM, limbal BM, and DM.*

Preliminary Data:

Preliminary in vitro data from our lab suggests that DM does indeed have some niche-like properties. Data from our labs shows that DM supports proliferation of LSCs in vitro. LSCs cultured on DM strongly maintain their stem cell phenotype (stemness) as well as their proliferative potential. LSCs cultured for 1 week on DM showed similar expression of putative LSC biomarkers p63 α (t-test, $p=0.34$) and ABCG2 (t-test, $p=0.29$) and similar co-localized expression of p63 α /ABCG2 (t-test, $p=0.20$) than cells grown on human amniotic membrane (HAM). This suggests that DM, as a substrate, is comparable to HAM at preserving the stemness of LSCs (Fig 1). Currently, HAM is the most widely used substrate for culturing and transplanting cultivated LSCs.²⁷ HAM is also widely used as a corneal on-lay alone in the treatment of ocular surface diseases such a persistent epithelial defects (PED) and LSCD.

On in-cell western, LSCs cultured for 1 week on DM showed comparable expression of putative LSC biomarkers: p63 α (t-test, $p=0.35$), ABCG2 (t-test,

$p=0.17$), and ABCB5 (t-test, $p=0.36$) when compared to HAM (Fig 2). And on BrdU pulse-chase assays, LSC cultured on DM showed similar proliferative rates to HAM (t-test, $p=0.52$) (Fig 3). This suggests that DM is comparable to HAM in maintaining both the stemness and proliferative potential of LSCs.




Figure 1. Limbal stem cells(LSC) isolated from donor tissue were seeded onto amniotic membrane (AM) and BrightMEM and cultured in serum-free media to confluence (10-14 days). LSC biomarker expression in cells cultured on BrightMEM was non-inferior and less variable than in cells cultured on AM.

Fig 2. Expression of LSC Biomarkers (In-cell Western)




Fig 2. Donor LSPC showed equal expression of LSC biomarkers (p63α, ABCG2, ABCB5) when cultured on DM versus HAM based on in-cell Western. ($\mu \pm SEM$).

Fig 3. BrdU Assay




Fig 3. BrdU assay showed similar proliferative rates in LSPC cultured on DM vs HAM ($\mu \pm SEM$).





Fig 4. Outgrowth culture of LSPCs on DM growing onto tissue culture plastic (TCP). Cells on DM (below the white arrows) show greater expression of p63α (red) and ABCG2 (green) compared to cells growing off DM on TCP.

Similarly, in outgrowth cultures of LSCs growing from DM onto adjacent tissue culture plastic (TCP), there was loss of ABCG2 and p63 α expression in cells that grew off DM after 2 weeks, while cells that remained supported on DM had greater retention of ABCG2 and p63 α expression (Fig 4). This supports the niche-like function of DM as a substrate.

LSCs cultured on DM also demonstrated normal stratification and differentiation into mature corneal epithelium (K12+) when airlifted, suggesting that DM can support normal stratified corneal epithelium as well as LSCs. Furthermore, in stratified epithelium cultured on DM, expression of ABCG2 and p63 α was preserved in the basal cell layer in contact with DM. This suggests that DM, as a substrate, can help preserve a reservoir of highly proliferative, LSC-like cells in the basal epithelium (Fig. 5).



In degradation studies, DM is also shown to be far more resistant to collagenase digestion than HAM. When submerged in serial dilutions of collagenase A (0.625 mg/ml, 1.25 mg/ml, 2.5 mg/ml, and 5 mg/ml), DM persisted for at least 24 hours at all concentrations, while HAM was completely degraded within 8 hours. Though HAM is effective in promoting corneal re-epithelialization, HAM is known to rapidly degrade within weeks on the ocular surface, limiting the durability of any therapeutic effect. Our preliminary data suggests that DM is far more durable and has the potential to be a long-term substrate on the ocular surface.

In early clinical cases, DM has been shown to adhere well and re-epithelialize in all cases (Fig. 6). DM has been shown to persist out to at least 13 months. In addition, early evidence suggests that DM may help prevent or delay recurrence of LSCD and corneal neovascularization in partial LSCD (Fig. 7).




Fig 6. POD#1 and POM#2 after DM corneal on-lay in two different patients. The backwards S-stamp orientation marker is seen at slit lamp confirming the presence of the graft under the epithelium.




Fig 7. Pre-op and POM#3 after DM corneal onlay + DSAEK for treatment of corneal endothelial decompensation and partial superior LSCD (2/2 trabeculectomy with MMC). Post-operatively, there was regression of the pre-operative superficial corneal neovascularization. Note: a bandage lens is present in the POM#3 image.

Summary:

Limbal niche components are potent drivers of LSC, TAC, and corneal epithelial proliferation and may have a therapeutic benefit in management of corneal epithelial disease. DM is naturally occurring basement membrane that is biochemically similar to limbal basement membrane, which is a critical component of the limbal niche. DM is routinely isolated for intraocular

transplantation by eye banks, optically clear, and resistant to collagenase digestion, making it an ideal long-term niche substitute on the ocular surface. Compared to HAM, which is widely used as a temporary corneal graft for treatment of ocular surface disease (HAM degrades after 1-2 weeks on the ocular surface), DM is relatively stable on the ocular surface and resistant to degradation. Preliminary data from our lab suggests that DM can provide a stable substrate on the corneal surface to support the proliferative potential of the basal epithelium. This suggests that DM may be an ideal corneal allograft for long-term support of LSCs and corneal epithelium on the corneal surface. Early clinical cases have demonstrated good adhesion and re-epithelialization of DM on the ocular surface. However, prospective clinical data is necessary for establishing the safety and efficacy of DM corneal allografts in patients with corneal epithelial disorders.

References:

1. Vemuganti, G.K., Sangwan, V.S. (2010). Interview: Affordability at the cutting edge: Stem cell therapy for ocular surface reconstruction. *Regen Medicine*, 5(3), 337-340.
2. Holland, E.J. (2015). Management of Limbal Stem Cell Deficiency: A Historical Perspective, Past, Present, and Future. *Cornea*, 34 (Supplement 10),S9-15.
3. Deng, S.X., Borderie, V., Chan, C.C., et al. (2019). Global Consensus on Definition, Classification, Diagnosis, and Staging of Limbal Stem Cell Deficiency. *Cornea*, 38(3), 364.
4. Chen, Z., de Paiva, C.S., Luo, L., Kretzer, F.L., Pflugfelder, S.C., Li, D.Q. (2004). Characterization of putative stem cell phenotype in human limbal epithelia. *Stem Cells*, 22(3), 355-366.
5. Zhou, M., Li, X., Lavker, R.M. (2006). Transcriptional profiling of enriched populations of stem cells versus transient amplifying cells. A comparison of limbal and corneal epithelial basal cells. *Journal of Biological Chemistry*, 281(28), 19600-19609.
6. Deng, S.X., Kruse, F., Gomes, J.A.P., et al. (2020). Global Consensus on the Management of Limbal Stem Cell Deficiency. *Cornea*, 39(10), 1291-1302.
7. Ballios, B.G., Weisbrod, M., Chan, C.C., et al. (2018). Systemic immunosuppression in limbal stem cell transplantation: best practices and future challenges. *Canadian Journal of Ophthalmology*, 53(4), 314-323.
8. Cheung, A.Y., Holland, E.J. (2017). Keratolimbal allograft. *Current Opinion in Ophthalmology*, 28(4), 377-381.
9. Sejpal, K., Yu, F., Aldave, A.J. (2011). The Boston keratoprosthesis in the management of corneal limbal stem cell deficiency. *Cornea*, 30(11), 1187-1194.
10. Koo, E.H., Hannush, S.B. (2021). Challenges in management of the Boston Keratoprosthesis Type 1. *Current Opinion in Ophthalmology*, 32(4).

11. Vaidyanathan, U., Hopping, G.C., Liu, H.Y., et al. (2019). Persistent corneal epithelial defects: A review article. *Medical Hypothesis Discovery and Innovation in Ophthalmology*, 8(3), 163-176.
12. Chan, E., Le, Q., Codriansky, A., Hong, J., Xu, J., Deng, S.X. (2016). Existence of Normal Limbal Epithelium in Eyes with Clinical Signs of Total Limbal Stem Cell Deficiency. *Cornea*, 35(11), 1483-1487.
13. Bonnet, C., Roberts, J.S., Deng, S.X. (2021). Limbal stem cell diseases. *Experimental Eye Research*, 205, 108437.
14. Huang, M., Wang, B., Wan, P., et al. (2015). Roles of limbal microvascular net and limbal stroma in regulating maintenance of limbal epithelial stem cells. *Cell and Tissue Research*, 359(2), 547-563.
15. Fukuda, K., Chikama, T., Nakamura, M., Nishida, T. (1999). Differential distribution of subchains of the basement membrane components type IV collagen and laminin among the amniotic membrane, cornea, and conjunctiva. *Cornea*, 18(1), 73-79.
16. Mariappan, I., Kacham, S., Purushotham, J., Maddileti, S., Siamwala, J., Sangwan, V.S. (2014). Spatial distribution of niche and stem cells in ex vivo human limbal cultures. *Stem Cells Translational Medicine*, 3(11), 1331-1341.
17. Lužnik, Z., Kopitar, A.N., Lapajne, L., et al. (2019). Identification and characterization of dendritic cell subtypes in preserved and cultured cadaveric human corneolimbal tissue on amniotic membrane. *Acta Ophthalmologica*, 97(2), e184-e193.
18. Mei, H., Gonzalez, S., Deng, S.X. (2012). Extracellular Matrix is an Important Component of Limbal Stem Cell Niche. *Journal of Functional Biomaterials*, 3(4), 879-894.
19. Kulkarni, B.B., Tighe, P.J., Mohammed, I., et al. (2010). Comparative transcriptional profiling of the limbal epithelial crypt demonstrates its putative stem cell niche characteristics. *BMC Genomics*, 11, 526.
20. Tseng, S.C.G., He, H., Zhang, S., Chen, S.Y. (2016). Niche Regulation of Limbal Epithelial Stem Cells: Relationship between Inflammation and Regeneration. *The Ocular Surface*, 14(2), 100-112.
21. Hirata, M., Yamaoka, T. (2018). Effect of stem cell niche elasticity/ECM protein on the self-beating cardiomyocyte differentiation of induced pluripotent stem (iPS) cells at different stages. *Acta Biomaterialia*, 65, 44-52.
22. Daya, S.M., Watson, A., Sharpe, J.R., et al. (2005). Outcomes and DNA analysis of ex vivo expanded stem cell allograft for ocular surface reconstruction. *Ophthalmology*, 112(3), 470-477.
23. Chen, P., Zhou, Q., Wang, J, et al. (2016). Characterization of the corneal surface in limbal stem cell deficiency and after transplantation of cultured allogeneic limbal epithelial cells. *Graefe's Archive for Clinical and Experimental Ophthalmology*, 254(9), 1765-1777.
24. Schlötzer-Schrehardt, U., Latta, L., Giebl, A., et al. (2021). Dysfunction of the limbal epithelial stem cell niche in aniridia-associated keratopathy. *The Ocular Surface*, 21, 160-173.
25. Ghiasi, M., Jadidi, K., Hashemi, M., et al. (2021). Application of mesenchymal stem cells in corneal regeneration. *Tissue and Cell*, 73, 101600.
26. Samaeekia, R., Rabiee, B., Putra, I., et al. (2018). Effect of human corneal mesenchymal stromal cell-derived exosomes on corneal epithelial wound healing. *Investigative Ophthalmology & Visual Science*, 59(12), 5194-5200.

27. Birbal, R.S., Sikder, S., Lie, J.T., Groeneveld-van Beek, E.A., Oellerich, S., Melles, G.R.J. (2018). Donor Tissue Preparation for Descemet Membrane Endothelial Keratoplasty: An Updated Review. *Cornea*, 37(1), 128-135.
28. Moolla, L., Mimouni, M., Din, N., et al. (2022). Effect of Collagenase A on Descemet Membrane Endothelial Keratoplasty Scroll Tightness. *Cornea*, 41(8), 1029.
29. van Dijk, K., Droutsas, K., Hou, J., Sangsari, S., Liarakos, V.S., Melles, G.R.J. (2014). Optical quality of the cornea after Descemet membrane endothelial keratoplasty. *Am Journal of Ophthalmology*. 158(1), 71-79.e1.
30. Dietrich-Ntoukas, T., Hofmann-Rummelt, C., Kruse, F.E., Schlötzer-Schrehardt, U. (2012). Comparative analysis of the basement membrane composition of the human limbus epithelium and amniotic membrane epithelium. *Cornea*, 31(5), 564-569.
31. Schlötzer-Schrehardt, U., Dietrich, T., Saito, K., et al. (2007). Characterization of extracellular matrix components in the limbal epithelial stem cell compartment. *Experimental Eye Research*, 85(6), 845-860.
32. Kabosova, A., Azar, D.T., Bannikov, G.A., et al. Compositional Differences between Infant and Adult Human Corneal Basement Membranes. *Investigative Ophthalmology & Visual Science*, 48(11), 4989-4999.
33. Echevarria, T.J., Chow, S., Watson, S., Wakefield, D., Di Girolamo, N. (2011). Vitronectin: a matrix support factor for human limbal epithelial progenitor cells. *Investigative Ophthalmology & Visual Science*, 52(11), 8138-8147.
34. Zhu, J., Wang, L.Y., Li, C.Y., et al. (2020). SPARC promotes self-renewal of limbal epithelial stem cells and ocular surface restoration through JNK and p38-MAPK signaling pathways. *Stem Cells*, 38(1), 134-145.